

SealMate™ SealPlate®

For fast and efficient microplate sealing

Application:

General laboratory use, short-term storage, tincubation and ELISA


Temperature Range:

-40°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
SM-KIT-SP	SealPlate Starter Kit, Dispenser + 2 Rolls	100/Roll
SM-KIT-SPS	SealPlate Starter Kit, Dispenser + 2 Rolls, Sterile	100/Roll
RSMSP-2	SealPlate Film, 2 Rolls	100/Roll
RSMSP-2-S	SealPlate Film, 2 Rolls, Sterile	100/Roll

The SealMate SealPlate Starter Kit features a re-usable dispenser and two rolls of perforated film. The dispenser facilitates the quick and easy dispensing of film from a roll allowing microplates to be sealed three times faster than using sheeted films. Cost-effective replacement rolls are sold separately. Ideal for use under the hood; the SealMate dispenser minimizes film contamination.

- Environmentally friendly; no paper liner to remove and discard; utilizes approximately 72% less material than SealPlate sheeted film
- Each sheet is separated by an adhesive-free, perforated, color-coded green band
- The smooth lower surface of the dispenser can be used to seal the film to the plate
- Dispenser protects unused seals from contamination
- Lower support prop prevents dispensing bill from contacting surfaces to prevent contamination
- Replacement rolls can be easily installed by releasing the locking clasps on the dispenser



"Check out SealPlate RollMate for 500 pre-scored sheets on a roll"